


Reliable, configurable, clear

Whether generated automatically or configured to meet your own requirements – the Vionex system continuously receives and processes all information on hygiene measures for the integrated Dentsply Sirona treatment centers. The hygiene operations are summarized in a clear workflow and analyzed intelligently. Vionex provides reliable support for complying with the strictest hygiene requirements in the dental clinic.



Dentsply Sirona
Sirona Dental Systems GmbH
Fabrikstraße 31, 64625 Bensheim, Germany
dentsplysirona.com

THE DENTAL
SOLUTIONS
COMPANY™

**Dentsply
Sirona**

Subject to technical changes and errors in the text. 0219

THE DENTAL
SOLUTIONS
COMPANY™

**Dentsply
Sirona**

Vionex

Unit management software for treatment centers
and simulation units

dentsplysirona.com/clinics



Procedural Solutions

- Preventive
- Restorative
- Orthodontics
- Endodontics
- Implants
- Prosthetics

Enabling Technologies

- CAD/CAM
- Imaging Systems
- Treatment Centers
- Instruments

Reliable hygiene management for every clinic

Vionex is a comprehensive software solution from Dentsply Sirona for dental clinics. It simplifies day-to-day hygiene and technical workflows and ensures a consistently high level of quality in hygiene documentation. It links all treatment centers and simulation workstations and saves data continuously and transparently on a central server. This allows the hygiene and technical data of each unit to be monitored at any time – quickly, clearly, completely – and from anywhere in the clinic.



Mobile comfort

Vionex users are mobile. Within the clinic network, every authorized user can work comfortably and quickly – check the hygiene status, document cleaning processes, or view service messages. Vionex can be accessed via a PC, tablet, or laptop – also via remote access.



Individual configuration

Vionex can be specifically oriented to different user groups. Customized configurations can be defined for each user group: dental assistants, hygiene officers, students, dentists, or service technicians. Special work methods, requirements and access rights can be taken into consideration. This results in personalized user interfaces; responsibilities are clearly assigned.



Intuitive operation

Simple and intuitive user interfaces and a clear menu navigation save time in workflows. The hygiene status and technical information of the treatment centers as well as customer specific floor plans are clearly visualized.

Your benefits with Vionex

Reliable infection prevention and optimal control

Vionex provides a quick overview and features useful tools for documenting and monitoring hygiene measures for treatment centers and simulation workstations.

Hygiene documentation

Purging and sanitizing processes carried out at the treatment center are documented in the software automatically and reliably. In addition, clinic-specific protocols can be defined, the manual implementation of which can be documented in detail. Step by step, this ensures a complete, digitally accessible overview of all hygiene measures implemented.

Hygiene reporting

The intelligent reporting function enables all data to be analyzed and presented in clear diagrams. Comprehensive filter and search functions can be used to highlight important parameters for conducting customized analysis. Added benefit: the analyses are a good basis for training the team in the correct implementation of hygiene guidelines.

Hygiene monitoring

The central data overview shows the status of all integrated treatment centers – completely and at a glance. Customized floor plans ensure better visual orientation. Proactive notifications by email if action is urgently required support the prompt implementation of hygiene measures.

Simulation Units

Vionex offers decisive advantages already in the skills lab:

- Vionex helps educating students of the right infection prevention measures already in pre-clinical semesters
- Vionex helps preventing errors by guiding students and newly-joined dental assistants through all hygiene measures to be carried out

Easy maintenance and fast service

Technical monitoring: fast, smart, and reliable

For service technicians, Vionex features an excellent overview of the current functioning of all units and tools for fast and easy servicing.



Benefits at a glance:

- access to service reports from all units at any time.
- allows you to be proactive: information displayed in the software enable early intervention so downtimes of treatment centers and simulation workstations can be minimized or even prevented. Preventive measures can be optimally planned and coordinated.
- assistance for technical servicing with an extensive library incl. instructions and solutions.
- comfortable remote access: within the clinic the technician can log on from any integrated PC, laptop, or tablet.