

Creating an environment with fewer potentially harmful aerosols



In a research study, **HVE removed 90% more aerosols** during ultrasonic scaling compared to a low-volume saliva ejector.⁴



DH20955 Purevac HVE Brochure 6pp A4.indd 5-1

Ordering Information

Item Number	Description
CD3006	Purevac* HVE Mirror Tips 6 Ultra FS

1. Al Maghlouth A, Al Yousef Y, Al Bagieh N. Qualitative and Quantitative Analysis of Bacterial Aerosols. J Contemp Dent Pract 2004 November; (5)4:091-100.

2. Milejczak CB. Optimum Travel Distance of Dental Aerosols in the Dental Hygiene Practice. J Dent Hyg. 2005 October;81(4):20-21

3. Christensen , R., PhD. (2006, November 1). Aerosols. Retrieved January 23, 2018, from http://www.dentistryiq.com/articles/wdj/print/volume-4/issue-10/you-and-your-practice/aerosols.htm

4. Jacks MJ: A laboratory comparison of evacuation devices on aerosol reduction. J Dent Hig 2002, 76, 202-206

For more information, visit dentsplysirona.com/purevac

@dentsplysirona.uk







Dentsply Sirona

Building 3, The Heights, Weybridge, Surrey KT13 ONY, United Kingdom

Telephone: +44 (0)1932 853 422

THE DENTAL SOLUTIONS COMPANY™













High Volume
Evacuation + Mirror
in a Single Instrument

dentsplysirona.com

Providing visibility and suction all in one hand

The new Purevac HVE Tip enables a one-handed approach to evacuating fluid and debris while facilitating retraction, visibility, and illumination during dental procedures.

Increases Efficiency

According to clinicians in a user study.*

Better Visibility


Fog-free mirror with continuous suction of fluids and debris.*

Greater Fluid Reduction

135% more than a low-volume saliva ejector.*

Reduced Noise

14% less than a standard HVE tip.* 5% less than a low-volume saliva ejector.*




Promoting better, safer, faster dental procedures

One innovative system that handles multiple procedures.









Polishing









Better visibility in new mirror:



Fog-free dental mirror provides indirect vision and illumination to the work site. Whilst the rounded and smooth edges are designed to aid retraction of the oral mucosa.

The Purevac HVE Mirror Tip is autoclavable for up to 100 cycles.



DH20955 Purevac HVE Brochure 6pp A4.indd 2-4