

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 1 of 7

1. Identification
Product identifier

MC Care Liquid

Relevant identified uses of the substance or mixture and uses advised against
Use of the substance/mixture

Cleaning agent

Details of the supplier of the safety data sheet

Company name: Sirona Dental Systems GmbH
 Street: Fabrikstraße 31
 Place: D-64625 Bensheim
 Telephone: +49 (0)625116-0
 e-mail (Contact person): <http://srvcontact.sirona.com/webformulars/EntryPage>
 Internet: www.dentsplysirona.com

Emergency telephone number: GBK (24 h) +49 (0)6132-84463

2. Hazard identification
Classification of the substance or mixture
WHMIS 2015

Serious eye damage/eye irritation: Eye Dam. 1

Label elements
WHMIS 2015
Signal word: Danger

Pictograms:

Hazard statements

Causes serious eye damage.

Precautionary statements

Wear protective gloves/protective clothing/eye protection/face protection.
 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.
 Continue rinsing.
 Immediately call a POISON CENTER/doctor.

Other hazards

No information available.

3. Composition/information on ingredients
Mixtures
Hazardous components

CAS No	Chemical name	Quantity
160875-66-1	Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-	5 - < 10 %
94441-92-6	Natrium-N-(2-carboxyethyl)-N-(2-ethylhexyl)-β-alaninat	1 - < 5 %
1554325-20-0	Quaternary C12-14 alkyl methyl amine ethoxylate methyl chloride	1 - < 5 %
64-17-5	ethanol, ethyl alcohol	1 - < 5 %

Safety Data Sheet

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 2 of 7

4. First-aid measures**Description of first aid measures****General information**

First aider: Pay attention to self-protection!

After inhalation

Provide fresh air. If unconscious place in recovery position and seek medical advice.

After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off contaminated clothing and wash it before reuse. In case of skin irritation, consult a physician.

After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water. Do NOT induce vomiting. Never give anything by mouth to an unconscious person or a person with cramps. Call a physician immediately.

Most important symptoms and effects, whether acute or delayed

Serious eye damage/eye irritation

Indication of immediate medical attention and special treatment needed

Treat symptomatically.

5. Fire-fighting measures**Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

Unsuitable extinguishing media

Full water jet

Specific hazards arising from the hazardous product

In case of fire may be liberated: Carbon dioxide (CO₂), Carbon monoxide.

Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Suppress gases/vapours/mists with water spray jet. Use water spray jet to protect personnel and to cool endangered containers. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

6. Accidental release measures**Personal precautions, protective equipment and emergency procedures**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

Environmental precautions

Do not allow to enter into surface water or drains. In case of gas escape or of entry into waterways, soil or drains, inform the responsible authorities.

Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

Reference to other sections

Safe handling: see section 7

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 3 of 7

Personal protection equipment: see section 8

Disposal: see section 13

7. Handling and storage
Precautions for safe handling
Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

Advice on protection against fire and explosion

No special fire protection measures are necessary.

Conditions for safe storage, including any incompatibilities
Requirements for storage rooms and vessels

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaustion at critical locations.

Advice on storage compatibility

Avoid: Acid, Strong alkali, Oxidising agent.

Further information on storage conditions

Keep away from heat. Protect from moisture.

8. Exposure controls/Personal protection
Control parameters
Exposure limits (ACGIH)

CAS No	Chemical name	ppm	mg/m ³	F/ml	Category	Origin
64-17-5	Ethanol	-	-		TWA (8 h)	ACGIH-2017
		1000			STEL (15 min)	ACGIH-2017

Exposure controls

Appropriate engineering controls

Provide adequate ventilation as well as local exhaustion at critical locations.

Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

Eye/face protection

Wear eye/face protection.

Hand protection

Suitable material:

HP Polyethylene (Breakthrough time (maximum wearing time): 10 - 60 min

NBR (Nitrile rubber) Thickness of the glove material: 0,4 mm

The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 4 of 7

Skin protection

Wear suitable protective clothing.

Respiratory protection

In case of inadequate ventilation wear respiratory protection.

Suitable material: Filter material/medium ABEK-P2.

Environmental exposure controls

Do not allow to enter into surface water or drains. In case of gas escape or of entry into waterways, soil or drains, inform the responsible authorities.

9. Physical and chemical properties**Information on basic physical and chemical properties**

Physical state:	Liquid
Colour:	greenish blue
Odour:	characteristic
pH-Value (at 20 °C):	6,6

Changes in the physical state

Melting point:	not determined
Initial boiling point and boiling range:	ca. 100 °C
Flash point:	not determined

Flammability

Solid:	not applicable
Gas:	not applicable

Explosive properties

No information available.

Lower explosive limits:	not determined
Upper explosive limits:	not determined
Ignition temperature:	not determined

Auto-ignition temperature

Solid:	not applicable
Gas:	not applicable

Decomposition temperature:	not determined
----------------------------	----------------

Oxidizing properties

Not oxidising.

Vapour pressure:	not determined
Density (at 20 °C):	1,015 g/cm ³
Water solubility:	easily soluble

Solubility in other solvents

not determined

Partition coefficient:	not determined
Viscosity / dynamic:	not determined
Viscosity / kinematic:	not determined
Vapour density:	not determined
Evaporation rate:	not determined

Other information

Odour threshold: No information available.

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 5 of 7

10. Stability and reactivity
Reactivity

No hazardous reaction when handled and stored according to provisions.

Chemical stability

The product is stable under storage at normal ambient temperatures.

Possibility of hazardous reactions

No known hazardous reactions.

Conditions to avoid

Humidity, Heat.

Incompatible materials

Avoid: Acid, Strong alkali, Oxidising agent.

Hazardous decomposition products

 In case of fire may be liberated: Carbon dioxide (CO₂), Carbon monoxide.

11. Toxicological information
Information on toxicological effects
Acute toxicity

Based on available data, the classification criteria are not met.

ATEmix tested

LD50, oral	Dose	Species	Source
	> 2000 mg/kg		calculated.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
160875-66-1	Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-				
	oral	ATE 500 mg/kg			
1554325-20-0	Quaternary C12-14 alkyl methyl amine ethoxylate methyl chloride				
	oral	LD50 > 300 -200 mg/kg	Rat	Manufacturer	

Irritation and corrosivity

Causes serious eye damage.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Sensitizing effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Name of toxicologically synergistic products

No information available.

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 6 of 7

12. Ecological information**Ecotoxicity**

The product is not: Ecotoxic.

Persistence and degradability

The product has not been tested.

Bioaccumulative potential

The product has not been tested.

Mobility in soil

The product has not been tested.

Other adverse effects

No information available.

Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

13. Disposal considerations**Waste treatment methods****Advice on disposal**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

14. Transport information**Canadian TDG****Proper shipping name:**

No dangerous good in sense of this transport regulation.

Marine transport (IMDG)**UN number:**

No dangerous good in sense of this transport regulation.

United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

Transport hazard class(es):

No dangerous good in sense of this transport regulation.

Packing group:

No dangerous good in sense of this transport regulation.

Air transport (ICAO-TI/IATA-DGR)**UN number:**

No dangerous good in sense of this transport regulation.

United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

Transport hazard class(es):

No dangerous good in sense of this transport regulation.

Packing group:

No dangerous good in sense of this transport regulation.

Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

15. Regulatory information**Canadian regulations****DSL/NDSL inventory status**

DSL:

according to WHMIS

MC Care Liquid

Revision date: 11.04.2018

Page 7 of 7

Quaternary C12-14 alkyl methyl amine ethoxylate methyl chloride: No
Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-: No
Natrium-N-(2-carboxyethyl)-N-(2-ethylhexyl)- β -alaninat: Yes.
ethanol, ethyl alcohol: Yes.
water: Yes.

NDSL:

Quaternary C12-14 alkyl methyl amine ethoxylate methyl chloride: No
Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-: No
Natrium-N-(2-carboxyethyl)-N-(2-ethylhexyl)- β -alaninat: No
ethanol, ethyl alcohol: Yes.
water: Yes.

National Pollutant Release Inventory (NPRI)

Quaternary C12-14 alkyl methyl amine ethoxylate methyl chloride: No
Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-: No
Natrium-N-(2-carboxyethyl)-N-(2-ethylhexyl)- β -alaninat: No
ethanol, ethyl alcohol: No
water: No

16. Other information**Abbreviations and acronyms**

DSL: Domestic Substances List
NDSL: Non Domestic Substances List
NPRI: National Pollutant Release Inventory
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonized System of Classification and Labelling of Chemicals
CAS: Chemical Abstracts Service
LC50: Lethal concentration, 50%
LD50: Lethal dose, 50%

Further Information

The information is based on present level of our knowledge. It does not, however, give assurances of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)