

+
**WE
KNOW
ENDO.**



X-Smart[®] Plus

Everything you like about X-Smart[®] with a plus

 **Dentsply
Sirona**
Endodontics

Everything you like about X-Smart® with a plus



Pre-programmed

For different file systems **Choose your preferred system at a single glance**

User friendly

Interface with «click and go» navigation and button controls dedicated to each function

No foot pedal



On/Off button on the motor handpiece

Excellent visibility and access

Due to the miniature contra-angle



Reciprocation

Single file root canal shaping with WaveOne® Gold the latest reciprocating technology to shape the canal with confidence.



Continuous motion

The next generation of the endo gold standard: ProTaper Next®, the performance of ProTaper® Universal with added versatility and a shorter clinical sequence.



The file library also includes preset programs for ProTaper® Gold, ProGlider®, PathFile®, Gates, Reciproc®, as well as 2 free programs for individual settings in continuous rotation.

Technical features

- Miniature contra-angle head that can be adjusted in 6 positions
- On/Off button on the motor handpiece
- Reciprocating motion and continuous rotation
- Includes a file library with preset programs for WaveOne®, WaveOne® Gold, ProTaper Next®, ProTaper® Gold, ProGlider®, PathFile®, Gates, Reciproc® – 2 free programs for individual settings in continuous rotation
- Auto reverse rotation at preset torque limit
- Torque range: 0.6 – 4.0 Ncm
- Speed range: 250 – 1'200 rpm
- Operates on a rechargeable battery and when connected to the mains
 - Battery autonomy: approx. 2 hours continuous use
 - Full recharge in approx. 5 hours

X-Smart® Plus



REF A103200000000

Endo motor with contra-angle and universal AC adaptor

X-Smart® Plus WaveOne® Gold Kit



REF B00XSPSW1GKIT

- X-Smart Plus Endo motor
- 5 blisters WaveOne Gold primary, 25mm (15 files)
- 4 blister WaveOne Gold assortment, 25mm (16 files)
- 1 blister WaveOne Gold Glider, 25mm (3 files)
- WaveOne Gold clinical procedure card

X-Smart® Plus 6:1 contra-angle



REF A103300000000

X-Smart® Plus cleaning adaptor



REF A100700000600

Cleaning and maintenance unit adaptor

X-Smart® Plus ProTaper Next® Kit



REF A103541000000

- X-Smart Plus Endo motor
- 10 blisters ProTaper Next assortment X1-X3, 25mm (30 files)
- 1 blister ProGlider, 25mm (3 files)
- ProTaper Next clinical procedure card