

# ProRinse® Endodontic Irrigation Probes

Effective Therapeutic Cleansing with safety on the side.




Our ProRinse probe is uniquely equipped with a lateral vent that gently stirs irrigants, while quickly and safely flushing them up and out of the canal. Unlike standard open-ended tips, ProRinse probe's lateral vents protect against extruding material periapically—a safety feature that is exclusive to the DENTSPLY line. In addition, because of its flexible stainless steel construction and 30-gauge tip, the probe easily delivers chemical cleaning solution at optimal depths.


Its distinct design delivers safe, effective therapeutic cleansing to the full length of the canal without risk to the apical foramen and facilitates instrument placement within the canal as filing progresses. The closed rounded tip and side-window port dispersal make ProRinse Endodontic Irrigation Probes ideal for use in canal irrigation.

Call your Tulsa Dental Specialties representative for more information on this exceptional product.

## THE PRO Rinse IRRIGATION PROBE ADVANTAGE:



- Vents solution laterally to ensure maximum irrigation and complete cleansing
- Delivers solution to optimal depth with 25mm length and 30-gauge tip
- Increases canal maneuverability through flexible design
- Offers excellent results with the highest degree of safety
- Features a Luer-connector for secure attachment and easy removal



See Directions For Use

For Dental Use Only

Rx only

**DENTSPLY Tulsa Dental Specialties**  
DENTSPLY International, Inc.  
5100 E. Skelly Drive, Suite 300  
Tulsa, OK 74135  
1-800-662-1202  
1-800-597-2779 (fax)  
[www.tulsadentalspecialties.com](http://www.tulsadentalspecialties.com)

**DENTSPLY**  
**TULSA DENTAL**  
**SPECIALTIES**  
Endodontics & Implants  
*Our Quality. Your Success.  
Their Satisfaction.*