

# Order form Primeprint Accessories and Solution Material



Company \_\_\_\_\_

Surname, First Name \_\_\_\_\_


Street, No. \_\_\_\_\_

Postal Code, City \_\_\_\_\_

Telephone \_\_\_\_\_ Fax \_\_\_\_\_

E-mail \_\_\_\_\_

Customer Number \_\_\_\_\_



## Accessories and consumables

Product	REF.	Quantity
Primeprint Activated Carbon Filter	6747823	
Collecting Container, Repl.	6757434	
Primeprint Box	6744895	
Primeprint Material Unit	6744903	
Primeprint Material Unit Type S <sup>1</sup>	6807965	
Inlay Transportation Container (5 pcs.)	6754241	
Foil with clamping frame spare	6744911	
S-Foil with clamping frame spare <sup>1</sup>	6807973	
Purosol optics cleaner	6782481	
Microfiber cloth	6802917	
Side Cutter DS	6758663	
Spatula DS	6758655	
Platform holder	6759402	
Washing Container 2.0	6819127	
Individual Exposure Carrier	6745074	
Primeprint PPU Activated Carbon Filter	6745603	
Primeprint PPU Ozone Filter	6752138	
Wiper Cards	6753631	

## Primeprint Solution Material

Product	REF.	Description	Medical device Class MDR	Medical device Class FDA	Quantity
Primeprint Splint	6740240	Resin for splints (high hardness compared to other splint resins for Primeprint); Filling weight: 1000 g	IIa	I	
Primeprint Splint SF <sup>2</sup>	6808856	Resin for splints, repositioners, mouth- & nightguards (high flexibility compared to other splint resins for Primeprint); Filling weight: 1000 g	IIa	I	
Primeprint Splint ST <sup>2</sup>	6808807	Resin for splints, repositioners, mouth- & nightguards (thinnest printable compared to other splint resins for Primeprint); Filling weight: 1000 g	IIa	I	
Primeprint Guide	6740232	Resin for surgical drilling guides; Filling weight: 1000 g	IIa	I	
Primeprint Temp A1	6740281	Resin for temporary crowns & bridges;; Filling weight: 1000 g	IIa	II	
Primeprint Temp A2	6740299				
Primeprint Temp A3	6740307				

<sup>1</sup> Required for printing of Primeprint Splint ST and Primeprint Splint SF.    <sup>2</sup> Only compatible with Material Unit Type S.

# Order form Primeprint Accessories and Solution Material



## Primeprint Solution Material

Product	REF.	Description	Medical device Class MDR	Medical device Class FDA	Quantity
Primeprint Model C	6824366	Resin for thermoforming, working and situation models; Filling weight: 1000 g	TEC resin	TEC resin	
Primeprint Model	6740257	Resin for working models; Filling weight: 1000 g	TEC resin	TEC resin	
Primeprint Model T	6740224	Resin for thermoforming models; Filling weight: 1000 g	TEC resin	TEC resin	
Primeprint Gingiva	6824291	Resin for dental gingiva masks; Filling weight: 500 g	TEC resin	TEC resin	
Primeprint Gingiva Starter Kit	6827609	Resin for dental gingiva masks bundled with two Washing Containers 2.0; Filling weight of the cartridge: 500 g	TEC resin	TEC resin	
Primeprint Cast	6740265	Resin for castables; Filling weight: 1000 g	TEC resin	TEC resin	
Primeprint Tray	6740216	Resin for impression trays; Filling weight: 1000 g	I	I	
Primeprint Lucitone Digital Print Original	0906401	Resin for denture bases; Filling weight: 1000 g	IIa	II	
Primeprint Lucitone Digital Print Light	0906402				
Primeprint Lucitone Digital Print Light Reddish Pink	0906403				
Primeprint Lucitone Digital Print Dark Reddish Pink	0906404				
Primeprint Lucitone Digital Value A1	0906405	Resin for try-ins; Filling weight: 1000 g	IIa	II	
Primeprint Lucitone Digital IPN 500ML A1	0906406	Resin for denture teeth; Filling weight: 500 g	IIa	II	
Primeprint Lucitone Digital IPN 500ML A2	0906407				
Primeprint Lucitone Digital IPN 500ML A3	0906408				
Primeprint Lucitone Digital IPN 500ML A3.5	0906409				
Primeprint Lucitone Digital IPN 500ML B1	0906410				
Primeprint Lucitone Digital IPN 500ML B2	0906411				
Primeprint Lucitone Digital IPN 500ML C1	0906412				
Primeprint Lucitone Digital IPN 500ML BL3	0906413				

## Denture Accessories

Product	REF.	Description	Quantity
LD Fuse Step 2 3D Bonding Resin OR 25g	0906186	Bonding resin to bond denture teeth to denture base	
LD Fuse Step 2 3D Bonding Resin LT 25g	0906188		
LD Fuse Step 2 3D Bonding Resin LRP 25g	0906189		
LD Fuse Step 2 3D Bonding Resin DRP 25g	0906190		
LD Fuse Step 3 Total 3D Sealer Size 15ml	0906399	Sealant that produces a smooth, glossy surface finish on the denture	
Biobond brush, #4 flat	9482646	Brush to apply bonding resin or sealer	