




DS Core[®]

Power your business growth with an integrated practice and efficient

- ✓ Enhance patient communication
- ✓ Expand service offering with new treatments
- ✓ Centralize collaboration with partners
- ✓ Connect equipment for efficient workflows
- ✓ Consolidate patient media and treatment files for efficient access



[Learn more](#)





Diagnose

Plan

Make

DS Core's wide range of features

<p>Coming soon in early 2025</p>	<p>Treatment Space Manage all treatment related files and workflows</p> <p>SureSmile® Simulator Illustrate the value of a SureSmile® aligner treatment</p> <p>CEREC Connected with DS Core Initiate CEREC design workflow and transfer data seamlessly¹</p>	<p>Manufacturing Manage all your manufacturing jobs in a single location²</p> <p>Ordering Manage all your lab and service orders in a single location</p>
<p>Online Media Library & Sharing Store and view patient media in one place, collaborate from anywhere</p>	<p>Canvas Enhance communication and engagement with your patients</p>	<p>DS Equipment Management Manage DS Core Powered and Connected equipment³</p>

Suscription overview

	Light <small>For small dental practices with an estimated usage of up to 90 DI scan and x-ray uploads / acquisitions per month.⁴</small>	Standard <small>For medium dental practices with an estimated usage of up to 270 DI scan and x-ray uploads / acquisitions per month.⁴</small>	Advanced <small>For large dental practices with an estimated usage of up to 540 DI scan and x-ray uploads / acquisitions per month.⁴</small>	Access <small>For practices that only want to send orders to labs.⁴ Connect Case Center replacement</small>
General Features				
Online Media Library	1 TB	6 TB	15 TB	Only temporary storage for up to 30 days ⁵
Sharing	✓	✓	✓	✗
Canvas	✓	✓	✓	✗
DS Equipment management with				
Automatic new IO-scan upload ⁶	✓	✓	✓	✓
Direct-to-cloud scanning with Primescan 2	✓	✓	✓	✗
Automatic new & historic x-ray upload ⁶	✓	✓	✓	Automatic new x-ray upload
3 rd party Equipment Support with				
Direct Upload from 3Shape Trios ⁷	✓	✓	✓	✗
Manual media upload	✓	✓	✓	✓
Plan Features				
Treatment Space	✓	✓	✓	✓
SureSmile® Simulator	✓	✓	✓	✗
CEREC connected with DS Core ⁸	✓	✓	✓	✗
Make Features				
Ordering	✓	✓	✓	✓
Manufacturing ⁹	✗	✓	✓	✗
EUR/month	19,95 € <small>Trial option available¹⁰</small>	99,95 €	229,95 €	At no cost

[→ Sign up on DSCore.com](https://www.ds-core.com)

¹ CEREC 5.3.1 license required; ² inLab CAM connected to DS Core required. ³ See "Dentsply Sirona Technology Solutions" section for an overview of "Powered by DS Core" and "Connected with DS Core" Equipment. ⁴ DS Core subscriptions support an unlimited number of users, however usage must align with the outlined recommendations from the order forms of DS Core per tier. ⁵ DS Core Access only provides temporary file storage. Uploaded files will be deleted after 7 days, if not used within an order. All files included in an order will be deleted 30 days once order is completed. ⁶ Automatic upload works with Connect/CEREC SW 5.2.3 or higher and Sidexis 4 V4.3.1 or higher. Use Connect 5.2.8, CEREC 5.3.1 and Sidexis 4 V4.4, for best performance and upload possibilities. ⁷ Seamless upload of IO scans from 3Shape Trios 3,4,5 scanners using 3shape Unite. ⁸ Works with CEREC 5.3.1 and higher. ⁹ Requires inLab CAM SW Connected with DS Core. ¹⁰ You can try all the features of DS Core Light for 3 months prior subscribing. You can terminate your DS Core Light subscription at any time during the first 3 months prior to your first payment being due. You will be asked to provide important feedback about your reasons for cancelling your subscription.

