

Scan bodies, model analogs and printed models – Compatibility chart

Atlantis FLO, Atlantis IO FLO and Elos Accurate IO Scan Body can be used with all connectivity when scanning for Atlantis abutments and are also included in the Atlantis Implant Kit for 3Shape (DME). 3Shape Scan Bodies and Elos Accurate Analog for Printed Models are also included in the DME file. This chart is valid for DME file created on 2019-10-01 or later.
















3Shape Scan Bodies, export the .stl files and send in the files utilizing Atlantis Scan Upload or any other .stl upload.

Atlantis Printed Model is available to order in Atlantis WebOrder and includes the Elos Accurate Model Analog for Printed Models.

Note: 3Shape Scan Bodies can only be used when sending in the scan files using the validated direct workflow from 3Shape to Atlantis WebOrder. The scan body detection must be done in the 3Shape software, therefore it is NOT possible to scan with



Manufacturer	Implant	Model scanning		Intraoral scanning			
		Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body	3Shape Scan Body <small>Only Atlantis Implant Kit for 3Shape (DME)</small>	Elos Accurate Analog for Printed Models <small>Only Atlantis Implant Kit for 3Shape (DME)</small>	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	REF	Supported
BioHorizons	Internal, Tapered Internal 3.5	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	BIY	PMA-ZSV35-1	YES
	Internal, Tapered Internal 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	BIG	PMA-ZSV45-1	YES
	Internal, Tapered Internal 5.7	B-03 (34529)	IO-B-03 (35214)	IO 5A-C	BIB	PMA-ZSV57-1	YES
	External Implant 3.5	A-03 (34518)	—	—	—	—	—
	External Implant 4.0	A-06 (34521)	—	—	—	—	—
	External Implant 5.0	A-09 (34524)	—	—	—	—	—
Biomet 3i	Certain 3.25, 4/3	F-01 (34539)	IO-F-01 (35225)	IO 7B-A	BCP	PMA-BCE34-1	YES
	Certain 4.0, 5/4	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	BCB	PMA-BCE41-1	YES
	Certain 5.0	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	BCY	PMA-BCE50-1	YES
	Certain 6.0	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	BCG	PMA-BCE60-1	YES
	Certain Prevail 3/4/3 (3.4), 4/3	F-01 (34539)	IO-F-01 (35225)	IO 7B-A	BCP	PMA-BCE34-1	YES
	Certain Prevail 4/5/4, 5/4	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	BCB	PMA-BCE41-1	YES
	Certain Prevail 5/6/5, 6/5	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	BCY	PMA-BCE50-1	YES
	Certain XP 4/5	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	—	PMA-BCE50-1	YES
	Certain XP 5/6	F-02 (34540)	IO-F-02 (35226)	IO 7B-B	—	PMA-BCE60-1	YES
	Miniplant 3.25	A-07 (34522)	—	—	BEP	—	—
	Osseotite 3.75, 4.0	A-07 (34522)	—	—	—	—	—
	Osseotite 5.0, 6.0	A-07 (34522)	—	—	—	—	—
	Osseotite XP 4/5	A-07 (34522)	—	—	—	BEY	—
	Osseotite XP 5/6	A-07 (34522)	—	—	—	BEG	—


Manufacturer	Implant	Model scanning		Intraoral scanning			
		Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body	3Shape Scan Body Only Atlantis Implant Kit for 3Shape (DME)	Elos Accurate Analog for Printed Models Only Atlantis Implant Kit for 3Shape (DME)	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	REF	Supported
Camlog	CONELOG 3.3	Q-01 (34666)	IO-Q-01 (35267)	—	—	—	—
	CONELOG 3.8	Q-02 (34667)	IO-Q-02 (35268)	—	—	—	—
	CONELOG 4.3	Q-02 (34667)	IO-Q-02 (35268)	—	—	—	—
	CONELOG 5.0	Q-03 (34668)	IO-Q-03 (35269)	—	—	—	—
	Screw-line Implant 3.3	N-01 (34584)	IO-N-01 (35239)	IO 9A-A	CCG	PMA-CAS33-1	YES
	Screw-line/Root-line Implant 3.8	N-02 (34585)	IO-N-02 (35240)	IO 9A-B	CCY	PMA-CAS38-1	YES
	Screw-line/Root-line Implant 4.3	N-03 (34586)	IO-N-03 (35241)	IO 9A-C	CCP	PMA-CAS43-1	YES
	Screw-line/Root-line Implant 5.0	N-04 (34587)	IO-N-04 (35242)	IO 9A-D	CCB	PMA-CAS50-1	YES
	Screw-line/Root-line Implant 6.0	N-04 (34587)	IO-N-04 (35242)	IO 9A-D	CCB	PMA-CAS60-1	YES
Dentsply Sirona <small>Note: Astra Tech Implant EV and OsseoSpeed EV have the same prosthetic interface.</small>	Ankylos C/X 3.5, 4.5, 5.5, 7.0 (non-indexed)	L-01 (34583)	IO-L-01 (35235)*	IO 8B-A SA	DYR Ø 3.5 DYY Ø 4.5 DYB Ø 5.5 DYG Ø 7.0	PMA-DAN35/70-1	YES
	Astra Tech Implant EV 3.0  GREEN	P-01 (34593)	IO-P-01 (35243)	IO 3B-A SA	DEG	PMA-ATE30-1	YES
	Astra Tech Implant EV EV 3.6  PURPLE	P-02 (34594)	IO-P-02 (35244)	IO 3B-B SA	DEP	PMA-ATE36-1	YES
	Astra Tech Implant EV EV 4.2  YELLOW	P-03 (34595)	IO-P-03 (35245)	IO 3B-C SA	DEY	PMA-ATE42-1	YES
	Astra Tech Implant EV EV 4.8  BLUE	P-04 (34596)	IO-P-04 (35246)	IO 3B-D SA	DEB	PMA-ATE48-1	YES
	Astra Tech Implant EV EV 5.4  BROWN	P-05 (34597)	IO-P-05 (35247)	IO 3B-E SA	DE5	PMA-ATE54-1	YES
	Astra Tech Implant EV Profile 4.2  YELLOW	P-06 (34598)	IO-P-06 (35248)	IO 3C-A SA	—	PMA-ATE42P-1	YES
	Astra Tech Implant EV Profile 4.8  BLUE	P-07 (34599)	IO-P-07 (35249)	IO 3C-B SA	—	PMA-ATE48P-1	YES
	OsseoSpeed / OsseoSpeed TX 3.0  YELLOW	C-01 (34530)	IO-C-01 (35215)	IO 3A-A SA	DAX	PMA-ATO30-1	YES
	OsseoSpeed / OsseoSpeed TX 3.5, 4.0  AQUA	C-02 (34531)	IO-C-02 (35216)	IO 3A-B SA	DAS	PMA-ATO35-1	YES
	OsseoSpeed / OsseoSpeed TX 4.5, 5.0  LILAC	C-03 (34532)	IO-C-03 (35217)	IO 3A-C SA	DAL	PMA-ATO45-1	YES
	OsseoSpeed TX Profile 4.5, 5.0	C-04 (34592)	IO-C-04 (35218)	—	—	—	—
	Xive S 3.0  030	J-01 (34575)	IO-J-01 (35230)	IO 8A-A	DX3	PMA-DXI30-1	YES
	Xive S 3.4  034	J-02 (34576)	IO-J-02 (35231)	IO 8A-B	DXG	PMA-DXI34-1	YES
	Xive S 3.8  038	J-03 (34577)	IO-J-03 (35232)	IO 8A-C	DXY	PMA-DXI38-1	YES
Xive S 4.5  045	J-03 (34577)	IO-J-04 (35264)	IO 8A-D	DXB	PMA-DXI45-1	YES	
Xive S 5.5  055	J-03 (34577)	IO-J-05 (35265)	IO 8A-D	DXR	PMA-DXI55-1	YES	
Hiossen	ET Mini 3.5	R-01 (34671)	IO-R-01 (35270)	—	OTM	—	—
	ET Regular 4.0, 4.5, 5.0, 6.0, 7.0	R-01 (34672)	IO-R-02 (35271)	—	OTR	—	—
Keystone Dental	PrimaConnex SD 3.3, Genesis 3.8	M-01 (34589)	—	—	—	—	—
	PrimaConnex SD 3.5	M-01 (34589)	—	—	—	—	—
	PrimaConnex RD 4.0, Genesis 4.5	M-02 (34590)	—	—	—	—	—
	PrimaConnex RD 4.1	M-02 (34590)	—	—	—	—	—
	PrimaConnex WD 5.0, Genesis 5.5, 6.5	M-03 (34591)	—	—	—	—	—








All trademarks and company names are the property of their respective owners.

* The Atlantis IO FLO for Ankylos C/X is for single use only.

** Astra Tech Implant System Lab Abutment Screw EV is designed for use with Astra Tech Implant System Implant Replica EV only and will not fit into the ELOS Accurate Model analog. Either order an extra Atlantis abutment screw, Abutment Screw EV or use ELOS lab screw.

*** Scanbody/Model analog exist but is not supported for Atlantis solutions.

Please turn over 

Manufacturer	Implant	Model scanning		Intraoral scanning			
		Atlantis FLO	Atlantis IO FLO	Elos Accurate IO Scan Body	3Shape Scan Body Only Atlantis Implant Kit for 3Shape (DME)	Elos Accurate Analog for Printed Models Only Atlantis Implant Kit for 3Shape (DME)	Atlantis Printed Model
		Code (REF)	Code (REF)	REF	REF	REF	Supported
MIS	M4 NP 3.3	B-04 (34678)	IO-B-04 (35273)	—	—	—	—
	M4 SP 3.75, 4.2	B-01 (34527)	IO-B-01 (35212)	—	—	—	—
	M4 WP 5.0, 6.0	B-02 (34528)	IO-B-02 (35213)	—	—	—	—
	SEVEN NP 3.3	B-04 (34678)	IO-B-04 (35273)	—	—	—	—
	SEVEN SP 3.75, 4.2	B-01 (34527)	IO-B-01 (35212)	—	—	—	—
	SEVEN WP 5.0, 6.0	B-02 (34528)	IO-B-02 (35213)	—	—	—	—
Nobel Biocare	NobelActive 3.0	K-03 (34665)	IO-K-03 (35266)	IO 2B-D SA	NA3	PMA-NBA30-1	YES
	NobelActive / NobelReplace Conical Connection NP 3.5 	K-01 (34578)	IO-K-01 (35233)	IO 2B-A SA	NAN	PMA-NBA35-1	YES
	NobelActive / NobelReplace Conical Connection RP 4.3 	K-02 (34579)	IO-K-02 (35234)	IO 2B-B SA	NAR	PMA-NBA43-1	YES
	NobelActive / NobelReplace Conical Connection RP 5.0 	K-02 (34579)	IO-K-02 (35234)	IO 2B-B SA	NAR	PMA-NBA43-1	YES
	NobelActive WP 5.5	K-04 (34677)	IO-K-04 (35272)	—***	NAW	PMA-NBA60-1	—
	NobelReplace NP 3.5  (Tri-channel connection)	D-01 (34533)	IO-D-01 (35219)	IO 2A-A	NRR	PMA-NBR35-1	YES
	NobelReplace RP 4.3  (Tri-channel connection)	D-02 (34534)	IO-D-02 (35220)	IO 2A-B	NRR	PMA-NBR43-1	YES
	NobelReplace WP 5.0  (Tri-channel connection)	D-03 (34535)	IO-D-03 (35221)	IO 2A-C	NRW	PMA-NBR50-1	YES
	NobelReplace 6.0  (Tri-channel connection)	D-04 (34536)	IO-D-04 (35222)	IO 2A-D	NR6	PMA-NBR60-1	—
	Brånemark System NP 3.3	A-04 (34519)	—	IO 6A-A	NBN	PMA-BRA35-1	YES
	Brånemark System RP 3.75, 4.0	A-06 (34521)	—	IO 6A-B	NBR	PMA-BRA41-1	YES
Brånemark System WP 5.0, 5.5	A-11 (34526)	—	IO 6A-C	NBW	PMA-BRA51-1	YES	
Straumann	Bone Level 3.3 NC (Narrow CrossFit)	H-01 (34542)	IO-H-01 (35228)	IO 4A-A SA	SBN	PMA-SBO33-1	YES
	Bone Level 4.1 RC (Regular CrossFit)	H-02 (34543)	IO-H-02 (35229)	IO 4A-B SA	SBR	PMA-SBO41-1	YES
	Bone Level 4.8 RC (Regular CrossFit)	H-02 (34543)	IO-H-02 (35229)	IO 4A-B SA	SBR	PMA-SBO41-1	YES
	Standard/Standard Plus 4.8 RN (Regular Neck)	E-01 (34537)	—	IO 4B-A	SSR	PMA-SSY48-1	—
	Standard/Standard Plus 4.8 WN (Wide Neck)	E-02 (34580)	—	IO 4B-B	SSW	PMA-SSY65-1	—
	Standard Plus 3.3 NN (Narrow Neck)	G-01 (34541)	—	—***	—***	—	—
Zimmer Dental	Screw-Vent / Screw-Vent friction fit 3.3, 3.7	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	—	PMA-ZSV35-1	YES
	Screw-Vent / Screw-Vent friction fit 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	—	PMA-ZSV45-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 3.5	B-01 (34527)	IO-B-01 (35212)	IO 5A-A	ZTG	PMA-ZSV35-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 4.5	B-02 (34528)	IO-B-02 (35213)	IO 5A-B	ZTP	PMA-ZSV45-1	YES
	Tapered Screw-Vent / Tapered Screw-Vent friction fit 5.7	B-03 (34529)	IO-B-03 (35214)	IO 5A-C	ZTY	PMA-ZSV57-1	YES
	SwissPluss 4.8	E-01 (34537)	—	—	—	—	—

All trademarks and company names are the property of their respective owners.

*** Scanbody/Model analog exist but is not supported for Atlantis solutions.