


Multiple-use drill usage tracking chart

Article	Number	Ø mm	Order no.		Usage									
					1	2	3	4	5	6	7	8	9	10
Guide Drill		-	68015133		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	①	1.9	68015324		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	②	2.35	68015325		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	③	2.95	68015326		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	④	3.55	68015327		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	⑤	4.15	68015328		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	⑥	4.75	68015329		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill short	⑦	5.35	68015330		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	①	1.9	68015331		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	②	2.35	68015332		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	③	2.95	68015333		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	④	3.55	68015334		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	⑤	4.15	68015335		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Drill long	⑥	4.75	68015336		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Tap		3.6	68015348		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Tap		4.2	68015349		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Tap		4.8	68015350		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Tap		5.4	68015351		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

Dispose of reusable drills and taps after a maximum of 10 reprocessing cycles/patients or if they are damaged or dull.

Multiple-use drill usage tracking chart

Article	Number	Ø mm	Order no.		Usage
Intermediate drills					1 2 3 4 5 6 7 8 9 10
Drill short	2½	2.65	68015337		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill short	3½	3.25	68015338		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill short	4½	3.85	68015339		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill short	5½	4.45	68015340		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill short	6½	5.05	68015341		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill long	2½	2.65	68015342		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill long	3½	3.25	68015343		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill long	4½	3.85	68015344		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill long	5½	4.45	68015345		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
Drill long	6½	5.05	68015346		☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐

Dispose of reusable drills and taps after a maximum of 10 reprocessing cycles/patients or if they are damaged or dull.