

Atlantis® suprastructures compatibility chart



Manufacturer/ Company	Implant system	Atlantis suprastructures prosthetic screw/ ASA screw availability	Screw Order REF number from Manufacturer/ Company	Torque Ncm
Dentsply Sirona	IMPLANT LEVEL			
	EV Connection: DS PrimeTaper Implant / DS OmniTaper Implant / Astra Tech Implant EV			
	EV (XS) Green*	x	25203	25
	EV (S) Purple	x	25204	25
	EV (M) Yellow	x	25205	25
	EV (L) Blue	x	25206	25
	EV (XL) Brown**	x	25207	25
	Astra Tech Implant EV Profile			
	Astra Tech Implant EV Profile 4.2	x	25205	25
	Astra Tech Implant EV Profile 4.8	x	25206	25
	Astra Tech Implant System			
	OsseoSpeed TX 3.0 - Yellow	x	24189	15
	OsseoSpeed TX 3.5, 4.0 - Aqua	x	24449	20
	OsseoSpeed TX 4.5, 5.0 - Lilac	x	24209	25
	OsseoSpeed TX Profile 4.5/5.0	x	24209	25
	OsseoSpeed 3.5 - Yellow (Previous)		24131	
	OsseoSpeed 4.0 - Blue (Previous)	x	24209	25
	Xive			
	Xive S 3.0	x	46-4304	24
	Xive S 3.4, 3.8, 4.5, 5.5	x	46-4305	24
	Frialit			
	Frialit 3.4, 3.8, 4.5, 5.5	x	46-4305	24
	IMZ			
	IMZ TwinPlus D3.3 Implant Level (After 1995)	x	54-1303	18
	IMZ TwinPlus D4 Implant Level (After 1995)	x	54-1303	18
	ABUTMENT LEVEL			
	EV Connection: DS PrimeTaper Implant / DS OmniTaper Implant / Astra Tech Implant EV			
	Multibase Abutment EV	x	26196	15
	UniAbutment EV	x	25481	15
	Angled Abutment EV	x (ASA not available)	25481	15
	Astra Tech Implant System			
	20° UniAbutment	x	22435	15
	45° UniAbutment	x	22435	15
	Angled Abutment 20°	x (ASA not available)	22435	15
	Ankylos			
	Ankylos Balance Base Narrow D4.2	x	3105 6022	10
	Ankylos Balance Base D5.5	x	3106 6022	10
	Xive			
	Xive TG 3.4, 3.8, 4.5	x	46-4300	24
	Xive MP 3.4, 3.8, 4.5	x	46-4300	24
	Xive MP 5.5	x	46-4305	24
	IMZ			
	IMZ TwinPlus D3.3 Abutment Level	x	51-1209	18
	IMZ TwinPlus D4 Abutment Level	x	51-1208	18

* Note: Not yet available for PrimeTaper EV

** Note: Not applicable for PrimeTaper EV

Please turn over →

Note: We recommend to check screw REF number and torque Ncm with the manufacturers, excluded Dentsply Sirona.

Note: Atlantis suprastructures assortment can vary from country to country. Please check Atlantis WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

Manufacturer/ Company	Implant system	Atlantis suprastructures prosthetic screw/ ASA screw availability	Screw Order REF number from Manufacturer/ Company	Torque Ncm
MIS	IMPLANT LEVEL			
	Internal Hex Narrow SEVEN - M4	x	MN-S0160	30
	Internal Hex Standard SEVEN - M4	x	MD-S0220	30
	Internal Hex Wide SEVEN - M4	x	MD-S0220	30
	Conical Connection Narrow C1 - V3	x	MN-S0160	30
	Conical Connection Standard C1 - V3	x	MD-S0220	30
	Conical Connection Wide C1	x	MD-S0220	30
	ABUTMENT LEVEL			
	Multi-Unit Abutment	x	MU-S0220	25
	CONNECT abutment D4.0	x	MM-S0160	30
	CONNECT abutment D5.7	x	MM-S0160	30
	Multiunit (Old)		MM-S0330	
BioHorizons	IMPLANT LEVEL			
	Internal/ Tapered 3.5, 4.5, 5.7	x	PXAS	30
	Internal/ Tapered/ Single 3.5, 4.5, 5.7	x	PXAS	30
	External Implant 3.5		130-300	
	External Implant 4.0, 5.0, 6.0		140-300	
	ABUTMENT LEVEL			
	Multi-unit abutment	x	PXMUPSR	15
	Internal D3.5/D4.5/D5.7 Abutment for screw		PXABS	
	D3.5/D4.5/D5.7 Angled Abutment for screw		PXABS	
	External D3.5/D4/D5 Abutment for screw		222-110	
Biomet 3i	IMPLANT LEVEL			
	Certain 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x	IUNIHT	20
	Certain Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	IUNIHT	20
	Certain XP 4/5, 5/6	x	IUNIHT	20
	MicroMiniplant 3.25	x	UNIHT	20
	Miniplant 3.25	x	UNIHT	20
	Osseotite 3.75, 4.0, 5.0, 6.0	x	UNIHT	20
	Osseotite XP 4/5, 5/6	x	UNIHT	20
	ABUTMENT LEVEL			
	Standard Abutment		GSH30	
	Conical Abutment D3.4		GSH30	
	Conical Abutment D4.1/D5		GSH30	
	TG Hex Abutment		GSH30	
	IOL Abutment		GSH30	
Low Profile Abutment	x	LPCTSH	10	
Camlog	IMPLANT LEVEL			
	Screw-line Implant 3.3	x	J4005.1601	20
	Screw-line/Root-line Implant 3.8, 4.3	x	J4005.1601	20
	Screw-line/Root-line Implant 5.0, 6.0	x	J4005.2001	20
	CONELOG D3.3, D3.8, D4.3 Implant Level		C4005.1601	
	CONELOG D5 Implant Level		C4005.2001	
	iSy Ø3.8/4.4/5.0		P4005.1614	
	ABUTMENT LEVEL			
	D3.3/D3.8/D4.3 Abutment level		J4005.1602	
	D5/D6 Abutment level		J4005.2002	
	COMFOUR Abutment System D3.3/D3.8/D4.3		J4012.1601	
	COMFOUR Abutment System D5/D6		J4012.2001	
	Vario SR Abutment D3.8/D4.3, D5/D6		J4005.2004	
	Vario SR Abutment D5/D6		J4005.2004	
Dyna	IMPLANT LEVEL			
	D3		82ES0	
	D3.6		31AS1	
	D4		31AS2	
	Octa Lock D4	x	81AS0	30
	Helix TM Octa	x	81AS0	30
	Octa Extension Level		82ES0	

Note: We recommend to check screw REF number and torque Ncm with the manufacturers, excluded Dentsply Sirona.

Please turn over →

Note: Atlantis suprastructures assortment can vary from country to country. Please check Atlantis WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

Manufacturer/ Company	Implant system	Atlantis suprastructures prosthetic screw/ ASA screw availability	Screw Order REF number from Manufacturer/ Company	Torque Ncm
Keystone Dental	IMPLANT LEVEL			
	PrimaConnex SD 3.3,3.5	x	45060K	30
	PrimaConnex RD 4.0,4.1	x	45060K	30
	PrimaConnex WD 5.0	x	45060K	30
	Genesis 3.8, 4.5, 5.5, 6.5	x	G71125	30
	Restore SD 3.3, 3.5		R9203-34-48K	
	Restore RD 3.75, 4.0		R9202-40-48K	
	Restore WD 5.0, 5.5, 6.0		R9203-50-48K	
	ABUTMENT LEVEL			
	MultiUnit RD		45164K	
	Stage1 RDS Abutment Level		S2408	
	Detachable Abutment		R9460-48K	
	Standard Abutment		R9200-40-48K	
Megagen	IMPLANT LEVEL			
	AnyOne 3.5	x	AS20	30
	AnyOne 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0	x	AS20	30
	ABUTMENT LEVEL			
	Multi-unit Abutment (N_Type)	x	MUAS	15
	Multi-unit Abutment (S_Type)	x	TASH140	15
	Octa Abutment 3.8	x	IRCS200	30
Octa Abutment 4.8	x	IRCS200	30	
Octa Abutment 5.8	x	IRCS200	30	
Neodent	IMPLANT LEVEL			
	GM Implant 3.5, 3.75, 4.0, 4.3	x	116.286	20
	GM Implant 5.0, 6.0, 7.0	x	116.286	20
	ABUTMENT LEVEL			
	GM Abutment	x	116266	10
GM Mini Conical Abutment	x	116269	10	
GM Micro Abutment	x	116269	10	
Nobel Biocare	IMPLANT LEVEL			
	NobelActive NP 3.5	x	36917	35
	NobelActive RP 4.3, 5.0	x	36918	35
	NobelActive WP 5.5	x	37892	35
	NobelReplace NP 3.5	x	29474	35
	NobelReplace RP 4.3	x	29475	35
	NobelReplace WP 5.0	x	29475	35
	NobelReplace 6.0	x	29475	35
	Branemark System NP 3.3	x	29282	35
	Branemark System RP 3.75, 4.0	x	29283	35
	Branemark System WP 5.0, 5.5		29284	
	ABUTMENT LEVEL			
	Multi-Unit Abutment RP	x	29285	15
	Multi-Unit Abutment WP		29286	
	Standard Abutment		32983	
	Estheticone Abutment		32981	
Steri-Oss PME Abutment		2346		
Hiossen/Osstem	IMPLANT LEVEL			
	GS/TS Mini		GSABSMT	
	GS/TS Standard		GSABSST	
	US Mini		ASM200	
	ABUTMENT LEVEL			
	GS Convertible Mini Abutment D4.0		GSFSM	
	GS Convertible Standard Abutment D4.8		GSFSR	
GS Convertible Wide Abutment D6.0		GSFSR		

Note: We recommend to check screw REF number and torque Ncm with the manufacturers, excluded Dentsply Sirona.

Please turn over 

Note: Atlantis suprastructures assortment can vary from country to country. Please check Atlantis WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

Manufacturer/ Company	Implant system	Atlantis suprastructures prosthetic screw/ ASA screw availability	Screw Order REF number from Manufacturer/ Company	Torque Ncm
Straumann	IMPLANT LEVEL			
	Bone Level 3.3 NC	x	025.2900	35
	Bone Level 4.1, 4.8 RC	x	025.4900	35
	NNC		048.314	
	Standard Plus 3.5 NN (Narrow Neck)		049.177	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x	048.356	35
	Standard/Standard Plus 4.8 WN (Wide Neck)	x	048.356	35
	BLX Implant RB 3.75, 4.5	x	065.0036	35
	BLX Implant WB 5.5, 6.5	x	065.0036	35
	ABUTMENT LEVEL			
	Monotype 8° Conus	x	048.350	15
	RN Abutment Level	x	048.350	15
	RN Abutment Level (Old)	x	048.350	15
	WN Abutment Level	x	048.350	15
	Bone Level Multi-Base Abutment D3.5, 4.5, 6.5	x	025.0900-04	15
	Bone Level Multi-Base Angled Abutment	x	025.0900-04	15
	Screw-Retained Abutment 3.5	x	023.4763	15
	Screw-Retained Abutment 4.6	x	023.4763	15
Zimmer Dental	IMPLANT LEVEL			
	Tapered Screw-Vent 3.5, 4.5, 5.7	x	MHLAS	30
	Screw Vent 3.3, 3.7, 4.5	x	MHLAS	30
	SwissPlus 3.8 Narrow	x	MHLAS	30
	SwissPlus 4.8		GPCAS	
	Advent D4.5		AVHLS	
	Spline 3.25, 4, 5mm		1537	
	ABUTMENT LEVEL			
	Screw Vent TSI Abutment		TSF	
	Tapered Abutment	x	SCTS	20
	Spectra Cone Abutment	x	SCTS	20
Spline Shouldered Abutment D4.5, 6.5		1621		

Note: We recommend to check screw REF number and torque Ncm with the manufacturers, excluded Dentsply Sirona.

Note: Atlantis suprastructures assortment can vary from country to country. Please check Atlantis WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.