

THE DENTAL
SOLUTIONS
COMPANY™



Sim Multi Bench

Training in dentistry and dental technology
in a compact space.



DENTISTRY

CEREC and prepCheck

Thanks to its flexible arm, an intraoral scanner can be used by several students during prosthodontics practice.



Ergonomic work
A healthy posture in every treatment position.



Work on animal specimens
Can also be used for exercises in surgery or periodontology.



Simply convert the phantom head into a surgery tray.

Training in dentistry and dental technology in a compact space

With Dentsply Sirona's customized bench solutions, you can give your students a realistic training environment enabling them to learn every discipline of dental technology. The Bench solutions can be extended with further optional dental add-ons to carry out interdisciplinary courses.

Compact

Make the most of the space you have thanks to its compact design.

Space-saving


Create space for as many students as possible.

Multi-functional

Conservative dentistry and prosthetic dental technology at a single workstation.

Always individually customized and matched to your needs.

Dentsply Sirona's simulation workstations are always planned with your input and built to suit your particular requirements and working conditions.



Use as a desk


Multimedia training system

Every student gets the best possible view and can keep up with the lesson.



Optimum use of space

Freedom of movement, even when carrying out dental technology tasks.



Overview of all options

LEDview Plus operating light

CEREC integration

Phantom head

Assistant element: high volume suction device and saliva ejector



prepCheck

Bunsen burner

Dry and wet grinding

Fully-fledged dentist element

Motor and handpiece for dental technology

The illustration shows a possible example of a fully equipped unit. The simulation workstations are always built to meet your own individual needs.

DENTAL TECHNOLOGY

Dry grinding

The protective disk and dry suction system prevent unhealthy dust from being inhaled.



Wet grinding

A basin built into the benchtop and a wet suction system allow practitioners to work without a spray mist.



Always a clean workplace

Thanks to the high volume suction device on the assistant element, the basin can always be cleaned quickly.



Dentsply Sirona

Sirona Dental Systems GmbH
Fabrikstraße 31, 64625 Bensheim, Deutschland
dentsplysirona.com

Procedural Solution

Preventive
Restorative
Orthodontics
Endodontics
Implants
Prosthetics

Enabling Technology

CAD/CAM
Imaging
Treatment Centers
Instruments