

Ceramco® 3 Veneering Porcelain

Inspire Your Passion

Ceramco 3

For more than 20 years, Ceramco 3 has been a proven high-fusing veneering system for metal-based frameworks, with a CTE of 13.9 to 15.1 $\mu\text{m}/\text{m} \cdot \text{K}$ (25-600 °C). Ceramco 3 allows dental technicians to have the freedom to work the way they work best.

- > Simplified layering system adaptable to various techniques and skill levels
- > Color-coded system provides a visual guide for build-up and shade accuracy
- > Controlled shrinkage results in confident layering and firing
- > Supports single to advanced build-ups and effects to inspire your passion



Inspire Your Passion with More

Opaque & Crystals

The base paste establishes the color foundation by masking the metal framework. The crystals help to create an intense light-scattering effect and improves the vitality of the restoration.

Classical & 3D-Shades


The Ceramco 3 Ceramic System is available in the traditional 16 A–D shades*, 26 Shade Series Dentines* keyed to the 3D-Master* shades and 8 illuminé bleach shades.

Build-up

The brilliant colorants assist in the build-up of restorations by reducing the training time and allowing technicians to clearly view the placement of each layer of porcelain.

Universal Stain & Glaze

For metal-based and all-ceramic restorations, eliminating multiple finishing systems in your lab.



Visit our website for more information at www.dentsplysirona.com/metalbasedcrownandbridge