



Classic®

## Mould Chart – Ultra Value Denture Teeth

For dental professionals and laboratories

---

THE DENTAL  
SOLUTIONS  
COMPANY™


 Dentsply  
Sirona



---

Everyone deserves the best denture they can afford. So, even for patients on a tight budget, a great looking denture is still possible. Classic Ultra Value denture teeth offer natural-like aesthetics, a wide selection of moulds and the most popular shades.

---



Classic®















# Anteriors Upper

Available in 24 shades and 1 bleach shade.

All dimensions in mm.  
Images scaled to size.









Mould No.  
(Lower  
Articulation)

Mould No.  
(Lower  
Articulation)





Mould No.  
(Lower  
Articulation)



Classic®









# Anteriors Lower

Available in 24 shades and 1 bleach shade.

All dimensions in mm.  
Images scaled to size.

Mould No.  
(Upper  
Articulation)

Mould No.  
(Upper  
Articulation)








Classic®





# Posterior Teeth

Available in 24 shades and 1 bleach shade.

## Classic 0° Posteriors








## Classic 10° Posteriors



**X** = The average width of a 1x4 set of upper posteriors  
**Y** = The average width of a 1x4 set of lower posteriors  
**U** = The average depth of the upper left first molar  
**L** = The average depth of the lower left first molar










ISO 22112: 2017. Processing instructions available upon request.

## Classic 20° Posteriors



All dimensions in mm.  
Images scaled to size.

## Classic 33° Posteriors



LS indicates Long Buccal-Short Lingual moulds.




# Classic® Anterior and Posterior Shades\*:

## 24 shades and 1 bleach shade

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, 59C, 62C, 65C, 66C, 67C, 69C, 77C, 81C, CW2.  
REF 905600

\*The A1-D4 designations correspond to the VITA classical A1-D4\* shade guide which is meant to be a guide, not a match. VITA classical A1-D4 is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.KG



## Denture Base

### Lucitone 199® Denture Base Resin



Lucitone 199 is the premium choice for both tissue and implant-supported full and partial dentures, offering outstanding aesthetics, high impact resistance, and superb flexural strength.

- Industry standard for premium high-impact denture base resins. Trusted and reliable.
- Superior high-impact resistance and flexural strength. Resists breakage.
- Balanced translucency and excellent vein simulation. Delivers aesthetic results.
- Compression packed, injected or microwavable if injected. Adapts to laboratory processes.
- Available in original, light pink, light reddish pink and dark pink shades. Covers patients' needs.

### Lucitone® HIPA High Impact Pour Acrylic

Lucitone HIPA is a high impact pourable acrylic denture base. It exceeds the ISO standards for impact resistance. Lucitone HIPA is a self-curing denture base material designed for use in fabrication, repair, rebasing, or relining of full and partial dentures including implant overdentures or other dental appliances. It is available in shades: Original, Light, Light Reddish Pink, Dark Pink, and Original Opaque.



Classic<sup>®</sup>  
Denture Teeth





# Combination Table

Select upper anterior mould form, and use table to determine recommended lower anterior and posterior options.

Classic® Anteriors		Classic® Posteriors			
Upper Mould No.	Articulates with Lower Mould	Combines with 33° Posterior Mould	Combines with 20° Posterior Mould	Combines with 10° Posterior Mould	Combines with 0° Posterior Mould
133	2C	30M	29M	F30	29M
134	3D	30M	31M	F33	31M
135	46	32M	33M	F33	33M
136	46	34M	33M	F33	33M
2D	2D	30L	29L	F32	29M
3D	3D	30M	31M	F33	31M
1N	46	30L	31M	F32	31M
2N	2N	30L	29L	F30	29M
2P	2P	30L	31L	F33	31M
A24	3N	30M	29M	F30	29M
A25	2E	30M	31M	F33	31M
A26	46	30M	31M	F33	31M
3M	3M	30M	29M	F30	29M
3N	3N	30L	31M	F32	31M
3P	3P	34M	33M	F33	33M
4M	3M	30M	29M	F30	29M
4N	2N	30L	31M	F30	31M
263	2N	30M	29L	F32	29M
264	2E	30M	31L	F33	31M
266	26	34L	33M	F33	33M
267	3R	34L	33M	F33	33M
268	3R	34L	31L	F33	33M

Dentsply Sirona | 800-243-1942 | [dentsplysirona.com](http://dentsplysirona.com)

© 2017 to 2021 Dentsply Sirona Inc. All rights reserved. Form No. 4033-A Rev. 18 (02/2021) 89564