


Sterilisers/Autoclaves

Sterilisation Pouches

Kanga Pouches

- Robust triple seal ensures strong seal without sacrificing the clean peel. Now with an added railing at the chevron for extra strength
- Multi-film construction improves film toughness and flexibility, prevents shattering
- Easy identification of contents through the blue tinted plastic film
- 60 gsm medical grade paper offers a superior barrier against contamination prior to use
- All imprints are positioned outside the actual packing area to prevent ink contamination reaching the product
- Quick and secure closing of a pack without a heat sealer
- Clean, fibre-free opening of packs ensures the safe presentation of sterile products
- Complete range of standard sizes for all surgery needs
- Suitable for use in hospital high vacuum and benchtop sterilisers



1057100 | Kanga Pouch 57 x 100 mm (pk 200)

1070230 | Kanga Pouch 70 x 230 mm (pk 200)

1090135 | Kanga Pouch 90 x 135 mm (pk 200)

1090230 | Kanga Pouch 90 x 230 mm (pk 200)

1135255 | Kanga Pouch 135 x 255 mm (pk 200)

Class B Autoclave Testing

Autoclave Consumables



AULT104 | Getinge Helix hollow load test set

Helix test for air removal and steam penetration in hollow loads (1 Helix Device & 250 Indicators)



AUPL254 | Getinge Bowie Dick Tests

Bowie Dick simulation test cards for air removal and steam penetration in porous loads (pack of 15 test cards)




AUS508 | Class 6 indicator test strips (pack of 250)

Getinge Titems Diagnostic GREEN Emulator

Replacement Trays and Racks

Standard Tray Rack



6486133 | DAC Professional Rack

For 5 open trays or 3 trays with lids

6532233 | DAC Professional Plus Rack

For 5 open trays or 3 trays with lids

6539790 | DAC Premium Rack

For 6 trays or 3 standard tray cassettes

6539808 | DAC Premium Plus Rack

For 6 trays or 3 standard tray cassettes

Nitra Rack for Instrument Packing



6117340 | Nitrarack 18L

For DAC Professional and DAC Premium improves the drying processed for stacking packaged instruments

6543222 | Nitrarack 22L

For DAC Professional Plus and DAC Premium Plus, improves the drying process of packaged instruments

Trays



6315290 | Small Autoclave Tray

29x19cm for sterilisers with a chamber depth of 35cm, recommended for DAC Professional and DAC Premium sterilisers.



MEAGTRAY | Large Autoclave Tray

42x19cm for sterilisers with a chamber depth of 45cm recommended for DAC Professional Plus and DAC Premium Plus



M011810 | MELASTore Tray 50

18 x 11.8 x 3 cm Sterilisable/Thermo Disinfectable Instrument Tray

M01181 | MELASTore Tray 100

27.5 x 17.6 x 3 cm Sterilisable/Thermo Disinfectable Instrument Tray

Printers and Data Storage

Patient Label Printing



668458 | MelaPrint 60

Barcode label printer for DAC Premium/DAC Premium Plus incl. a roll of labels

6543180 | Labels Nitraprint 60, white

Adhesive labels for barcodes

1 roll with 1,000 labels, double-sided adhesive

6543198 | Labels Nitraprint 60, blue

Adhesive labels for barcodes

1 roll with 1,000 labels, double-sided adhesive

Data Log Printing



GDOTPRINTER | Dot Matrix Printer

Printer Compatible with DAC Professional and DAC Universal

GPROSIRSCABLE | Printer Cable

Printer Cable for DAC Professional

GDACCABLE I Printer Cable

Printer Cable for DAC Universal (S)

GPLUGPACK | Printer Power Supply



6051770 | Dac Universal S Printer

6599018 | Paper for DAC Universal S (pk 5)



M01144 | MELAprint 44 Dot Matrix Printer

Printer compatible with DAC Universal & DAC Professional. Can also be utilized with DAC Premium if combined with M40295 Ethernet Module

Datalogging



6543214 | Nitraflash CF Card

For data collection on all devices

9145155 | Lexar USB 3.0 Card Reader

For the movement of data on DAC Premium and DAC Premium Plus devices

6134394 | Nitraflash Data Logger

Direct connection to DAC Professional

JMP-3144 | 5V DC Adaptors for Nitraflash

For the storage of Data on DAC Professional Devices

DAC Universal

DAC Universal Handpiece Oil



6259118 | NitramOil #2 Bottle: blue, 200 mL
Maintenance concentrate for the DAC Universal with following serial numbers:
14.000-99.999
104.000-199.999
202.000-999.999
Also suitable for DAC Universal S

DAC Universal Cleaning Tablets



6635499 | NitraClean cleaning tablets
For the cleaning of the chamber (50 pcs)

Filters



6078583 | DAC Universal filter (6 pcs)
6263045 | DAC Universal Replacement Sterile filter
6504489 | DAC Universal S filter



6698166 | DAC Universal S Drain Filter (6x)
6742519 | O-Ring for Drain Filter small (10x)
6742501 | O-Ring for Drain Filter large (10x)

DAC Universal Testing Process Challenge Device



DAC Universal PCD Test Body **6051820**
DAC Universal S PCD Test Body **6742956**

DAC Universal Indicators (also suitable for other Class S Sterilisers)



3M1243B | 3M comply Class 5 Sterigage strips
Steam Chemical integrators for pack control monitoring for all 118-138°C (245-280°F) steam sterilisation cycles.
•Recommended for DAC Universal and other Class S Sterilisers



6742857 | DAC Universal S Chemical Indicator Strips (510x)

DAC Universal Indicator Holders

6051788 | DAC Universal Indicator Holder (standard lid)
6542489 | DAC Universal Indicator Holder (flex lid)
6743624 | DAC Universal S Indicator Holder

DAC Universal S Lids



6742907 | DAC Universal S Pink Lid
For the sterilisation of dental handpieces, without adaptors.



6742931 | DAC Universal S White Lid
For the sterilisation of ultrasonic instruments and multifunctional syringes, without adaptors.



6709815 | DAC Universal Blue Lid
For the disinfection only of dental handpieces, without adaptors.



6709823 | DAC Universal S Green Lid

For the disinfection only of ultrasonic instruments and multifunctional syringes, without adaptors.



6709997 | DAC Universal S Check and Clean Lid



6710003 | DAC Universal S Check and Clean Cap

DAC Universal Lids



6051739 | DAC Universal Standard Lid

For the sterilisation of dental handpieces, without adaptors.



6535855 | DAC Universal Flex Lid

For the sterilisation of ultrasonic instruments and multifunctional syringes, without adaptors.

6085612 | Gasket Screws for Adapters Standard (6x)

6535921 | O-Ring screws for Adapters Flex (6x)

DAC Universal Lid Holders



N720024 | DAC Universal Lid Holder

6709856 | DAC Universal S Lid Holder

Water Treatment Accessories

MELAdem 47 Filters



M37470 | Mixed Bed Resin Filter MELAdem47

M37460 | Carbon Filter MELAdem47

M37450 | Fine Filter MELAdem

MELAdem 40 Filter

M61026 | Filter MELAdem 40

Nitradem Filter

66932098 | Nitradem Mini Filter

ANZ ✉ clientservices@dentsplysirona.com 📍 11-21 Gilby Road, Mount Waverley, VIC 3149 🌐 www.dentsplysirona.com

Australia ☎ 1300 552 929 | **New Zealand** ☎ 0800 336 877

YOUR
CONFIDENCE
IS OUR
INSPIRATION.

**Dentsply
Sirona**