

IPN 3D™

Portrait Inspired

Digital Denture Teeth



Mould Cross References

Anterior Upper

IPN 3D Anterior Upper	Portrait IPN Mould Inspirations	ATS Facial Form Inspiration
110DP	13E, 21D	Square, Square Tapering
112HP	21D, 21C	Square Tapering
120JP	32E, 22E	Square Ovoid, Square Tapering
122MP	12E, 32E	Square, Square Ovoid
142QP	21X, 22G	Square, Square Tapering
163SP	12G, 32G	Square, Square Ovoid
185VP	11H, 21G	Square, Square Tapering

Anterior Lower

IPN 3D Anterior Lower	Portrait IPN Mould Inspirations
5DP	H
5HP	F, C
5JP	E, H
5MP	G, N
5QP	P, L
5SP	R, V, I
5VP	S, Rx, K1

Posterior

IPN 3D Posterior Upper	IPN 3D Posterior Lower	IPN 3D Size	Portrait IPN Mould Inspirations
7L10P	9L10P	L	334
7M10P	9M10P	M	332
7S10P	9S10P	S	330
7L33P	9L33P	L	34M
7M33P	9M33P	M	32M
7S33P	9S33P	S	30M

IPN 3D and Portrait IPN mould forms are not exact matches.