

Midwest® E Plus

Electric Handpiece Attachments

Available
September 2016

PRECISION & CONTROL AT YOUR FINGERTIPS

Perform a broad range of procedures with the power and flexibility provided by the Midwest E Plus Electric Handpiece Attachments. They easily let you apply your skills to the fullest with confidence, accuracy, and efficiency.

- Smaller head eases access*
- Titanium shell reduces weight*
- Ceramic head bearings extend durability*

MORE POWER. MORE TORQUE. LESS NOISE. LESS WEIGHT.

Driving the Midwest E Plus Handpiece Attachments to their full capabilities is the dynamic Midwest E motor, weighing just 59 grams and measuring only 1.25". It's one of the smallest, lightest electric motors you'll find today, yet it delivers tremendous, consistent power and torque.

With this motor, you can efficiently handle even the toughest crown and bridge procedures, and manage multi-unit cases faster. Plus, it's whisper-quiet, creating a more relaxing, soothing environment for patients and staff alike. Smartly designed, the Midwest E motor includes:

- Up to 3.0Ncm of torque
- LED lighting for increased visibility
- Adjustable light intensity and light delay through color LCD touchscreen
- A speed range from 100-40,000 RPM
- Short-nosed E-Type connection for wide compatibility

*As compared to original Midwest E attachments.
Ceramic head bearings in 1:5 attachment only.



ONE-TOUCH CONTROL

Using the intuitive LCD color touchscreen controls, you can quickly view and select the proper settings for high-speed, low-speed, and endodontic procedures. Customizing the attachment operation from one screen means you can complete procedures easily and efficiently.

Along with adjusting the speed, torque limits, and light operation, the touchscreen:

- Has four user-defined presets in both “Prep” and “Endo” functions
- Is easily customizable with 28 pre-programmed gear ratios and 2 user-defined ratios
- Allows for seamless integration into existing dental units with detachable faceplate and remote mounting kit
- Endo functionality includes user-defined torque limits and 3 auto-reverse options



KEEP IT CLEAN

In clinical studies, correct cleaning and lubrication have proven to be the most important factors for handpiece performance and longevity. To ensure you get the most out of your handpieces, the Midwest® Automate® system:

- Reduces handling and cleanup time by more than 80%
- Automatically removes debris, precisely lubricates all moving parts, and accurately expels excess
- Ensures consistent maintenance every time
- Ends the excess mess of manual handpiece maintenance
- Optimizes long-term performance with quality maintenance



COMPLETE CUTTING SOLUTIONS

As a leader in both handpiece and bur manufacturing, only Midwest offers the ultimate integration of design, performance, and integrity.



ORDER INFORMATION

Midwest® E Electric Handpiece System

876000 Midwest® E Electric Handpiece System

Midwest® E Plus Attachments

875205 Midwest® E Plus 1:5 High Speed Attachment

875210 Midwest® E Plus 1:1 Low Speed Attachment

875220 Midwest® E Plus 6:1 Endo Attachment

875230 Midwest® E Plus 1:1 Straight Attachment

System Components

875080 Midwest® E Remote Mounting Kit

Midwest® Automate®

Automated Maintenance Unit

380200 Midwest® Automate® Unit

380520 Absorption Pad (6 Pack)

380080 Midwest Plus™ Aerosol Spray Can

Warranties

Midwest E Electric Handpiece System is backed by nearly 90 years of Midwest handpiece experience, including:

3-year warranty on motor, tubing and control unit

2-year warranty on attachments



Customer Service: 1-800-800-2888

901 West Oakton Street
Des Plaines, IL 60018-1843 USA

www.dentsplysirona.com