

Midwest[®] E Pro

Electric Handpiece Attachments

3-Year
Warranty
Protection

Smaller Size. Longer Life.

Perform a broad range of procedures with the power and flexibility provided by the Midwest E Pro Electric Handpiece Attachments.

- Smaller and lighter design improves ergonomics*
- Diamond-like carbon-coated bearings increase longevity*
- 3-year warranty reduces repair costs



More Power. More Torque. Less Noise. Less Weight.

Drive the Midwest E Pro Handpiece Attachments to their full capabilities with the dynamic Midwest E motor. Weighing just 59 grams and measuring only 1.25", it's one of the smallest, lightest electric motors you'll find today.

With this motor, you can efficiently handle even the toughest crown and bridge procedures, and manage multi-unit cases faster. Plus, it's whisper-quiet, creating a more relaxing, soothing environment for patients and staff alike. Smartly designed, the Midwest E motor includes:

- Up to 3.0Ncm of torque
- LED lighting with adjustable intensity
- Short-nosed E-Type connection for wide compatibility

*As compared to original Midwest E attachments.



One-Touch Control

Using the intuitive LCD color touchscreen controls, you can quickly view and select the proper settings for high-speed, low-speed, and endodontic procedures. Customizing the attachment operation from one screen means you can complete procedures easily and efficiently.

Along with adjusting the speed, torque limits, and light operation, the touchscreen:

- Has four user-defined presets in both “Prep” and “Endo” functions
- Is easily customizable with 28 pre-programmed gear ratios and 2 user-defined ratios
- Allows for seamless integration into existing dental units with detachable faceplate and remote mounting kit
- Endo functionality includes user-defined torque limits and 3 auto-reverse options



Keep it Clean

In clinical studies, correct cleaning and lubrication have proven to be the most important factors for handpiece performance and longevity. To ensure you get the most out of your handpieces, the Midwest® Automate® system:

- Reduces handling and cleanup time by more than 80%
- Automatically removes debris, precisely lubricates all moving parts, and accurately expels excess
- Ensures consistent maintenance every time
- Ends the excess mess of manual handpiece maintenance
- Optimizes long-term performance with quality maintenance



Complete Cutting Solutions

As a leader in both handpiece and bur manufacturing, only Midwest offers the ultimate integration of design, performance, and integrity. Explore our full array of carbide and diamond burs with the [Midwest Bur Advisor](#), a virtual bur wallet and utilities program available for Apple, Android and Windows. Download it free at the Apple Store, on Google Play, or through www.MidwestBurApp.com.



Order Information

Midwest® E Pro Attachments *3-year warranty protection*

- 875105 Midwest® E Pro 1:5 High Speed Attachment
- 875110 Midwest® E Pro 1:1 Low Speed Attachment

Midwest® E Electric Handpiece System

Motor, tubing and control unit *3-year warranty protection*

- 876000 Midwest® E Electric Handpiece System

Midwest® E System Components

- 875080 Midwest® E Remote Mounting Kit

Midwest® E Plus Attachments *2-year warranty protection*

- 875205 Midwest® E Plus 1:5 High Speed Attachment
- 875210 Midwest® E Plus 1:1 Low Speed Attachment
- 875220 Midwest® E Plus 6:1 Endo Attachment
- 875230 Midwest® E Plus 1:1 Straight Attachment

Midwest® Automate® Automated Maintenance Unit

- 380200 Midwest® Automate® Unit
- 380520 Absorption Pad (6 Pack)
- 380080 Midwest Plus™ Aerosol Spray Can