

Your individual options

Here you can view the complete list of individual options, and a checklist of what is included with basic units.

Intego Pro

Patient chair		Pro Ambidextrous
2 treatment positions, S and O position	■	■
Integrated junction box	■	■
Chair lift system 185 kg	■	■
Flat headrest	■	■
Double articulating headrest	□	□
Premium upholstery	■	■
Lounge upholstery	□	□
ErgoMotion	□	□
4-way foot switch	■	■
Booster cushion for small patients	□	□
Right arm rest	□	□
Left arm rest	□	□
Valve set for external device	□	□
Vacuum pressure reduction kit	□	□
Cover installation kit	□	□
Separate junction box	□	□
Adapter plate	□	□

Dentist element		Pro Ambidextrous
Height-adjustable support arm with pneumatic brake	■	■
5 instrument positions	■	■
Additional camera holder for TS dentist element	□	—
Flexible TS camera holder	—	□
Hydro-pneumatic supply for 2 drill drives w/o fiber optics	■	■
Pneumatic instrument foot control	■	■
C+ electronic foot control	□ S, A, I	□ S, A, I
EasyPad user interface with 2 user-specific settings	■	■
EasyTouch user interface with 4 user-specific settings	□ A, I	□ A, I
Sprayvit E with light	□ S, A, I	□ S, A, I
Sprayvit E with light (water button on the left)	□ S, A, I	□ S, A, I
Turbine hose with light	□ S, A, I	□ S, A, I
1st BL ISO E motor with hose	□ S, A, I	□ S, A, I
2nd BL ISO E motor with hose	□ I	□ I
SiroSonic L scaler	□ A, I	□ A, I
Third drill drive	□ I	□ I
Satelec Mini LED curing light	□	□
Intra-oral light source kit for turbine and Sprayvit E	■	■
Torque control for endodontic treatments	□ E	□ E
eControl: Speed range of 100-40.000 rpm	□ E	□ E
Integrated ApexLocator	□ E	□ E
Endo 6:1 contra-angle handpiece	□ E	□ E
Intra-oral x-ray viewer	□	□
Purge and Auto Purge	■	■
Flexible TS tray holder for 2 standard size trays	—	□
Additional TS tray holder for 2 standard size trays	□	—
Silicone pad for TS dentist element	□	□
Additional CS tray holder for 1 standard size tray	□	□
Additional CS tray holder for 2 standard size trays	□	□
Silicon pad for CS dentist element	□	□
Removable silicone sleeves for handles	□	□

Assistant element		Pro Ambidextrous
Compact assistant element, swivel mounted directly on the cuspidor arm	—	■
Comfort assistant element, ideal for two-handed dentistry	■	□
4 instrument positions	■	■
Saliva ejector and HVE	■	■
Sprayvit E with light	□ I	□ I
Sprayvit E with light (water button on the left)	□ I	□ I
2nd HVE kit	□	□
Satelec Mini LED curing light	□	□
Touch pad	□	□
Light source for Sprayvit E	■	■
Silicone pad for Comfort assistant element	□	□

Water unit		Pro Ambidextrous
Prepared for integrated disinfection device	■	■
Cuspidor with manual swivel function	■	■
White cuspidor	■	■
Glass cuspidor	□	—
Water unit without cuspidor	□	□
Multi-station valve	■	■
Cuspidor valve for standard wet suction	□	□
Automatic separator	□	□
Amalgam separator	□	□
Integrated disinfection system with permanent water-line disinfection and automatic sanitization function	□	□
Water heater	■	■
Water flush suction hose cleaning system	■	■
Chemical suction hose cleaning system	□	□

Operating light/Heliodont Plus options		Pro Ambidextrous
LEDlight	□ S	□ S
LEDview Plus	□ I, T	□ I, T
LEDview Plus ceiling light	□	—
Adaptation for LEDview Plus/Heliodont Plus ceiling combination	□	—
Adaptation for Heliodont Plus unit model	□	—

Patient communication		Pro Ambidextrous
USB interface in the dentist or assistant element	□ SiR, SiS	□ SiR, SiS
SiroCam F	□ SiS	□ SiS
SiroCam AF	□	□
SiroCam AF+	□	□
Monitor mount on the light support column	□ SiS	□ SiS
22" AC monitor (with integrated power supply unit)	□ SiS*	□ SiS*
Network interface	□ SiR, SiS	□ SiR, SiS

■ Included in basic unit
 □ Option
 — Not possible/applicable
 SiR Sivation Ready package
 SiS Sivation Start package
 SiS* Sivation Start package with monitor

S Start Pro package
 A Advance Pro package
 I Innovation Pro package
 T Touch Plus package
 E Endo Pro package

Availability of options and features can vary by country.