


Atlantis® abutments implant compatibility chart



Manufacturer	Implant	Atlantis Healing Abutment		Atlantis Abutment			Atlantis Crown ²		Atlantis CustomBase solution - Atlantis Abutment & Atlantis Crown ²		Atlantis Crown Abutment		Atlantis Conus Abutment	
							Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Abutment material												
		Titanium Gold-shaded ti	Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia ¹	Titanium Gold-shaded ti	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium		
Dentsply Sirona	INTERNAL													
	EV Connection													
	PrimeTaper EV 3.6 - (S) Purple	x	x		x		x	x	x	x		x	x	
	PrimeTaper EV 4.2 - (M) Yellow	x	x		x		x	x	x	x		x	x	
	PrimeTaper EV 4.8 - (L) Blue	x	x		x		x	x	x	x		x	x	
	PrimeTaper EV 5.4 - (L) Blue	x	x		x		x	x	x	x		x	x	
	OmniTaper EV 3.0 - (XS) Green	x	x		x		x	x	x	x				
	OmniTaper EV 3.4 - (S) Purple	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	OmniTaper EV 3.8 - (M) Yellow	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	OmniTaper EV 4.5 - (L) Blue	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	OmniTaper EV 5.5 - (XL) Brown	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	Astra Tech Implant EV 3.0 - (XS) Green	x	x		x		x	x	x					
	Astra Tech Implant EV 3.6 - (S) Purple	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	Astra Tech Implant EV 4.2 - (M) Yellow	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	Astra Tech Implant EV 4.8 - (L) Blue	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	Astra Tech Implant EV 5.4 - (XL) Brown	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	Astra Tech Implant EV Profile													
	Astra Tech Implant EV Profile 4.2 Yellow	x	x	x	x		x	x	x			x	x	
	Astra Tech Implant EV Profile 4.8 Blue	x	x	x	x		x	x	x			x	x	
	Astra Tech Implant System TX													
	OsseoSpeed/OsseoSpeed TX 3.0 - Yellow	x	x		x		x	x	x					
	OsseoSpeed/OsseoSpeed TX 3.5, 4.0 - Aqua	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	OsseoSpeed/OsseoSpeed TX 4.5, 5.0 - Lilac	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	OsseoSpeed TX Profile 4.5, 5.0	x	x		x		x		x			x	x	
	OsseoSpeed 4.0 - Blue (Previous)		x										x	
	Xive / Frialit													
	Xive S 3.0	x	x		x		x	x	x					
	Xive S 3.4, 3.8, 4.5, 5.5	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
Frialit 3.4, 3.8, 4.5, 5.5	x	x	x	x	x ¹	x	x	x	x	x	x	x	x	
Ankylos														
Ankylos C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)		x	x	x								x		
Ankylos C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x	x	x	x ¹	x	x				x	x	x	
MIS	MIS													
	INTERNAL													
	C1 NP 3.3	x	x		x		x	x	x			x	x	
	C1 SP 3.75, 4.2	x	x		x		x	x	x			x	x	
	C1 WP 5.0	x	x		x		x	x	x			x	x	
	M4 NP 3.3	x	x		x		x	x	x			x	x	
	M4 SP 3.75, 4.2	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	M4 WP 5.0, 6.0	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	SEVEN NP 3.3	x	x		x		x	x	x			x	x	
	SEVEN SP 3.75, 4.2	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	SEVEN WP 5.0, 6.0	x	x	x	x	x ¹	x	x	x	x	x	x	x	x
	V3 NP 3.3	x	x		x		x	x	x			x	x	
	V3 SP 3.9, 4.3, 5.0	x	x		x		x	x	x			x	x	

All trademarks and company names are the property of their respective owners.

Please turn over →

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.

Manufacturer	Implant	Atlantis Healing Abutment		Atlantis Abutment		Atlantis Crown ²		Atlantis CustomBase solution – Atlantis Abutment & Atlantis Crown ²		Atlantis Crown Abutment		Atlantis Conus Abutment	
						Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Abutment material											
		Titanium Gold-shaded ti	Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia ¹	Titanium Gold-shaded ti	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium	
BioHorizons	BioHorizons												
	INTERNAL												
	Internal 3.5, 4.5, 5.7	x	x	x	x	x ¹	x	x	x	x	x	x	x
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	x ¹	x	x	x	x	x	x	x
	EXTERNAL												
	External implant 3.5		x ¹										x ¹
External implant 4.0, 5.0		x ¹	x							x		x ¹	
Biomet 3i	Biomet 3i												
	INTERNAL												
	Certain 4/3, 5/4	x	x	x	x	x ¹	x			x	x	x	x
	Certain 3.25, 4.0, 5.0, 6.0	x	x	x	x	x ¹	x			x	x	x	x
	Certain Prevail 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	x ¹	x			x	x	x	x
	Certain XP 4/5, 5/6	x	x	x	x	x ¹	x			x	x	x	x
	EXTERNAL												
	Miniplant 3.25		x	x						x	x	x	x
Osseotite 3.75, 4.0, 5.0, 6.0		x	x						x	x	x	x	
Osseotite XP 4/5, 5/6		x	x						x	x	x	x	
Camlog	Conelog												
	INTERNAL												
	CONELOG 3.3	x	x		x		x			x		x	x
	CONELOG 3.8, 4.3	x	x		x		x			x		x	x
	CONELOG 5.0	x	x		x		x			x		x	x
	Screw-line Implant 3.3	x	x		x		x			x		x	x
Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x ¹	x			x	x	x	x	
Keystone Dental	Keystone Dental												
	INTERNAL												
	PrimaConnex SD 3.3, 3.5, Genesis 3.8	x	x	x	x	x ¹	x			x	x	x	x
	PrimaConnex RD 4.0, 4.1, Genesis 4.5	x	x	x	x	x ¹	x			x	x	x	x
PrimaConnex WD 5.0, Genesis 5.5, 6.5	x	x	x	x	x ¹	x			x	x	x	x	
Nobel Biocare	Nobel Biocare												
	INTERNAL												
	NobelActive 3.0	x	x		x		x	x	x				
	NobelActive/NobelReplace Conical Connection NP 3.5	x	x	x	x	x ¹	x	x	x	x	x	x	x
	NobelActive/NobelReplace Conical Connection RP 4.3, 5.0	x	x	x	x	x ¹	x	x	x	x	x	x	x
	NobelActive WP 5.5	x	x	x	x	x ¹	x	x	x	x	x	x	x
	NobelReplace NP 3.5 (Tri-channel connection)	x	x	x	x	x ¹	x			x	x	x	x
	NobelReplace RP 4.3 (Tri-channel connection)	x	x	x	x	x ¹	x			x	x	x	x
	NobelReplace WP 5.0 (Tri-channel connection)	x	x	x	x	x ¹	x			x	x	x	x
	NobelReplace 6.0 (Tri-channel connection)	x	x	x	x	x ¹	x			x	x	x	x
	EXTERNAL												
Brånemark System NP 3.3		x		x		x			x		x	x	
Brånemark System RP 3.75, 4.0		x	x	x		x			x	x	x	x	
Brånemark System WP 5.0, 5.5		x	x	x		x			x	x	x		
Osstem	Osstem												
	INTERNAL												
	Osstem TS Mini 3.5	x	x	x	x		x			x	x	x	x
Osstem TS Regular 4.0, 4.5, 5.0, 6.0, 7.0	x	x	x	x		x			x	x	x	x	

All trademarks and company names are the property of their respective owners.

Please turn over 

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.

Manufacturer	Implant	Atlantis Healing Abutment	Atlantis Abutment		Atlantis Crown ²		Atlantis CustomBase solution – Atlantis Abutment & Atlantis Crown ²		Atlantis Crown Abutment		Atlantis Conus Abutment	
					Cut-back or Full-contour		Cut-back or Full-contour		Cut-back		Custom	Over-denture
		Abutment material										
		Titanium Gold-shaded ti	Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia ¹	Titanium Gold-shaded ti	Angulated screw access	Titanium	Zirconia	Titanium Gold-shaded ti	Titanium
Straumann	Straumann											
	INTERNAL											
	Bone Level 3.3 NC (Narrow CrossFit)	x	x	x	x	x ¹	x	x	x	x	x	x
	Bone Level 4.1, 4.8 RC (Regular CrossFit)	x	x	x	x	x ¹	x	x	x	x	x	x
	Standard/Standard Plus 4.8 RN (Regular Neck)	x	x		x		x		x		x	x
	Standard/Standard Plus 4.8 WN (Wide Neck)	x	x		x		x		x		x	x
	EXTERNAL											
Standard Plus 3.3 NN (Narrow Neck)		x ¹								x ¹		
Zimmer Dental	Zimmer											
	INTERNAL											
	Tapered Screw-Vent 3.5, 4.5, 5.7	x	x	x	x	x ¹	x	x	x	x	x	x
	Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	x	x		x		x		x		x	x
	Screw-Vent 3.3, 3.7, 4.5	x	x	x	x	x ¹	x	x	x	x	x	x
	Screw-Vent 3.3, 3.7, 4.5 friction fit	x	x		x		x		x		x	x
	SwissPlus 4.8	x	x		x		x		x		x	x

All trademarks and company names are the property of their respective owners.

¹ Manufactured in Waltham, US. Turnaroundtime: 5-10 extra working days (excluding US Shipping).

² Includes zirconia crowns and Atlantis Temporary Crowns in PMMA.